










Common causes of Dry Eye include:

-  • Sun and wind exposure
-  • Cold or dry air
-  • Indoor heating and air conditioning
-  • Prolonged computer screen viewing
-  • High altitudes
-  • Eye surgery
-  • Hormonal fluctuations


Many medications can contribute to Dry Eye including:

-  • Antihistamines and decongestants
- Blood pressure medication
- Anti-depressants and anti-anxiety medications

Additionally many pre-existing conditions or diseases also effect the eye's ability to lubricate, for example:

- Thyroid deficiencies
- Sjogren's Syndrome
-  • Rheumatoid Arthritis
- Autoimmune disorders
- Bell's Palsy
- Diabetes

While severity of conditions vary from mild, moderate to extreme, the common symptoms include:


-  • Redness
- Burning
- Itching
- Scratchiness
- Excessive tearing
- Foreign body in the eye


WE ARE HAPPY WHEN YOU EXPERIENCE DRY EYE RELIEF LIKE NEVER BEFORE.

“ I never felt compelled to write a company to thank them for a product – but your product is so wonderful and has improved the quality of my life, that I had to let you know.” – **LISA, CA**

“ Hurray for Tranquileyes®! Using your product has helped protect my eyes from the normal progression of damage due to Sjogren's. You have a lifelong customer. Keep up the good work.” – **SUSAN, FL**

“ Last night I wore my Tranquileyes® for the first time. WOWOWOWOW!!! Amazing! How wonderful it felt to wake up and NOT have my eyes dry and achy. Thank you! – **JOANNE, IL** ”

 **EYE ECO™ GIVES BACK.** Since day one, we have donated a percentage of our sales to the Sjogren's Syndrome Foundation in an effort to help find a cure for this serious autoimmune disease, of which dry eye is a hallmark symptom. This is who we are and will always be.

 To learn more about Sjogren's, visit www.sjogrens.org



ADDRESS 28544 Old Town Front Street, Suite 303
Temecula, California 92590 USA
PHONE 951.695.4393 or 1.888.730.7999
FAX 951.695.3393
WEB eyeeco.com



NATURAL SOOTHING SOLUTIONS FOR IMMEDIATE AND LONG-TERM DRY EYE RELIEF.

ABOUT OCULAR SURFACE WELLNESS

The eye's inability to produce sufficient tears and lubricate itself is a condition that affects over 40 million people every day.

The causes of Dry Eye Syndrome include aging, environmental irritants, lifestyle, injuries or medical conditions.

Treating dryness is not only important for your comfort, but also for your eye's health. Besides lubricating the eye, the tear film helps fight infection, provides nourishment, and creates a smooth surface on the cornea, keeping your vision clear.

Though some people suffer from all of the symptoms associated with Dry Eye Syndrome, many experience only a few. If you experience any of the symptoms, ask your doctor about Dry Eye. There is currently no cure for Dry Eye Syndrome, but there are a variety of treatments to address the symptoms including the application of moist-heat. For more than ten years, Eye Eco™ has been the leading provider of moist-heat options for Dry Eye relief.




GOOD BETTER BEST

MOIST-HEAT/COLD THERAPIES

STUDIES HAVE PROVEN THAT APPLYING MOIST-HEAT TO THE EYELIDS RELIEVES DRY EYE.

Temperature and duration are important to effectively loosen clogged oils which slow tear evaporation, relieving dryness. Equally important is the humidity. Increased humidity stabilizes the tear film, promotes circulation and healing while, hydrating sensitive eyelid skin. Eye Eco™ offers three scientifically proven solutions for dry eye relief. All are safe, natural, re-usable and easy to clean. Using a hygrometer to measure the temperature and relative humidity, these tests indicate the Tranquileyes® system maintains the ideal balance for optimum performance.

SIMPLY PUT, THE LONGER THE OPTIMUM LEVELS OF MOIST-HEAT, THE BETTER.



WHICH eyeeco SOLUTION IS BEST FOR YOU?


- derm™ DRY EYE RELIEF MASK
- tranquileyes® EYE HYDRATING THERAPY
- tranquileyes® xl EYE HYDRATING THERAPY
- thermoeyes™ Instant
- thermoeyes™ Instant XL
- thermoeyes™ Beads
- thermoeyes™ Beads XL

For more information, ask your doctor or visit eyeeco.com.




The Dry Eye Relief Mask, D.E.R.M.™, is our entry level solution for immediate dry eye relief. The mask contains Silica beads which naturally attract moisture from the air and release it when microwaved. Heating the mask for 20 seconds will provide soothing moist-heat relief for 4-6 minutes. Applying a moistened cotton liner to the mask prior to heating will extend the treatment to 7-9 minutes. D.E.R.M.™ is hypoallergenic and can also be used to create cold therapies for allergy relief.

derm™ 5-6 minutes moist-heat



derm™ 7-9 minutes moist-heat with Moist Liner




Our patented Tranquileyes® system, combines the goggle, warm water and reusable gel packs to create soothing moist-heat therapies which sustain the optimum temperature, duration and humidity required for immediate and long-term relief. Also comfortable for nighttime shelter.


Tranquileyes® and Tranquileyes® XL systems can be used to create 15 minute cold therapies, helping to relieve allergies, puffiness and migraines.



thermoeyes™ Instant 18-20 minutes moist-heat




thermoeyes™ Beads 12-15 minutes moist-heat



Our Tranquileyes® XL goggle is thirty percent larger than the original tranquileyes. The larger area accommodates a larger gel pack to sustain optimum moist-heat therapies up to 25 minutes along with extended coverage. The Tranquileyes® XL is recommended for immediate and long-term relief from severe dry eye, MGD (Meibomian Gland Dysfunction), Blepharitis and ocular allergies.

thermoeyes™ Instant XL 20-25 minutes moist-heat



thermoeyes™ Beads XL 20-25 minutes moist-heat

